



The first and only
single-use, 14-day*, continuously recording ECG monitor

The Zio® Patch is indicated for use on patients who may be asymptomatic or who may suffer from transient symptoms such as:

- **Palpitations**
- **Anxiety**
- **Syncope**
- **Light-headedness**
- **Dizziness**
- **Fatigue**
- **Pre-syncope**
- **Shortness of breath**


FEATURES

- Single-use
- Up to 14 days of ECG recording*
- Small, lightweight, water-resistant
- 100% of every device is recycled

**Total wear time may vary by patient*

BENEFITS

- The extended monitoring time provided by the Zio® Patch beyond the 24 – 48 hours of monitoring provided by traditional Holter monitors may significantly improve diagnostic yield
- Small size can increase patient compliance
- Single-use devices may limit the risk of transmitting infections between patients



DEVICE SPECIFICATIONS

TECHNICAL FEATURES

Simple 1-button operation to mark symptoms

MEMORY CHARACTERISTICS

ECG Channels	1 channel
Memory capacity	14 days
Recording Format	Continuous

ELECTRICAL CHARACTERISTICS

Frequency Response	0.15 Hz to 34 Hz
Input Impedance	> 3 Mohm
Differential Range	+/- 1.65 mV
Resolution	10 bits

PHYSICAL CHARACTERISTICS


Dimensions	123 x 53 x 10.7 mm
Weight	34 grams
Shelf Life	6 months

POWER REQUIREMENTS

Battery Type	Lithium Polymer
Battery Life	14 days

STANDARDS COMPLIANCE

AAMI/ANSI EC38
IEC 60601-1



iRhythm Technologies, Inc., the provider of Zio® products and services, markets Zio® in collaboration with St. Jude Medical, Inc.

